

CLIP V CIDE—DISINFECTANT, FUNGICIDE, GERMICIDE & VIRUCIDE

Product Description:

CLIP-V-CIDE is for use in nursing homes, health care facilities, tanning facilities, pedicure foot spa units, beauty salons, barbershops and hairstyling establishments.

Features:

Use on beauty/barber instruments and tools, hard non-porous surfaces such as floors, sinks, telephones, and toilet seats. All raw materials used for production undergo a thorough quality control process before being used for manufacturing high quality products.

Specifications:

Color:	Amber
Odor:	None Added
pH:	2
Physical State:	Liquid

Available In:

5 Gallon Pail	Item# 98PXT05
30 Gallon Pails	Item # 98PXT30
55 Gallon Pails	Item# 98PXT55
4x1 Gallon Bottles	Item# 98PXT0G

Directions:

SURFACES: On inanimate, hard non-porous surfaces such as floors, sink tops, tables, acrylic tanning beds, telephones, toilet seats, and door handles. Dilute 2 ounces of CLIP-V-CIDE in 32 ounces of water. Apply prepared solution to pre-cleaned surfaces with a cloth, mop, sponge, or coarse mist pump or trigger sprayer so as to wet all surfaces thoroughly. Allow to remain wet for 10 minutes, then let air-dry. **BEAUTY/BARBER INSTRUMENTS AND TOOLS:** Thoroughly clean brushes, combs, clipper blades, shears, razors, and other inanimate hard surfaced non-porous tools and implements prior to complete immersion for 10 minutes (or as required by local authorities) in a solution as prepared above. Rinse to remove CLIP-V-CIDE after immersion and dry before use. Fresh solution must be prepared daily or more often if solution becomes diluted or soiled. **PERFORMANCE:** When tested with full immersion for 10 minutes in a solution as prepared above, CLIP-V-CIDE has been found to be effective against: *Pseudomonas aeruginosa*, *Salmonella enterica*, *Staphylococcus aureus* (including *Staphylococcus aureus* Methicillin Resistant (MTSA)), *Trichophyton mentagrophytes*, *Herpes Simplex Virus Type-1 and Type-2, *Human Immunodeficiency Virus Type 1 (HIV-1)(AIDS virus), Influenza A Virus (in 5% Serum), *Hepatitis B Virus and *Hepatitis C Virus. **SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATING AGAINST HIV-1, HBV OR HCV ON SURFACES/OBJECTS SOILED WITH BLOOD AND BODY FLUIDS.** On pre-cleaned environmental surfaces/objects previously soiled with blood/body fluids in health care settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood or body fluids and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 (HIV-1) (associated with AIDS), Hepatitis B Virus (HBV), or Hepatitis C Virus (HCV).

First Aid:

Call poison control center or doctor immediately for treatment advice. **IF IN EYES:** Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing eye. **IF ON SKIN OR CLOTHING:** Take off contaminated clothing. Rinse skin immediately with plenty of water for 15-20 minutes. **IF SWALLOWED:** Have person sip a glass of water if able to swallow. **DO NOT** induce vomiting unless told to do so by the poison control center or doctor. Do not give anything by mouth to an unconscious person. Have the product container or label with you when calling a Poison Control Center or doctor or going for treatment. **NOTE TO PHYSICIAN:** Probable mucosal damage may contraindicate the use of gastric lavage.

Storage:

Store in original container protected from direct sunlight in a dry, cool and well-ventilated area, away from incompatible materials and food and drink. Use appropriate containment to avoid environmental contamination. Keep container tightly closed and sealed until ready to use. Containers that have been opened must be carefully resealed and kept upright to prevent leakage. Do not store in unlabeled containers. Use appropriate containment to avoid environmental contamination.



**INNOVATIVE
CHEMICAL
CORPORATION**

7769 95th St. S, Cottage Grove, MN 55016
P:651.649.1762 TF:866.472.3932 F:651.649.1787
www.iccmn.com